

# ELASTIC TRUSSES AND OTHER PRODUCTS



Code - 601

Code - 602

- Used before and after operations for hernia.
- Prevents the development and appearance of hernia before operations and the loosing of the stitches after operations.
- It is comfortable, thanks to its special designed pad.



nia. ce ng

**Elastic Truss (Double)** 

**Silver Socks** 

Code - 602

• Right and left types are available.



- Used for all types of sports.
- Avoids the hanging of the testis, gives support those which hang.
- Also used wounded testis, hemorrhage and varicose or infection of testis.



NEW

Recommended especially for • diabetic patients.

Inner part is made from silver yarn. •

# Thanks to the natural specialty of the silver:

- Speeds the treatment up by accelerating blood circulation.
  Prevents the bad bacterium to generate, and removes the odor.
  - Available in two sizes. Foot Nr: 36-40 Small Size

Foot Nr: 41-45 Large Size





Code - 975

Therapeutic



ELLY SA

Silver Tubular Bandage

Code - 400



### Advantages of the silver yarn:

- Speeds the treatment up by accelerating blood circulation.
- Prevents the bad bacterium to generate, and removes the odor.
  - It is used at the limbs like hand, wrist, foot, leg, arm.





# **ELASTIC TRUSSES AND OTHER PRODUCTS**



### **Elastic Bandages**

• Used in order to support arm, wrist, foot, leg and other places.

Code - 406 Elastic Bandage 6 cm x 2 m

Code - 408 Elastic Bandage 8 cm x 2 m

Code - 410 Elastic Bandage 10 cm x 2 m

Code - 412 Elastic Bandage 12 cm x 2 m

Code - 415 Elastic Bandage 15 cm x 2 m

## Code - 871

### **Sport Short**

- Maximizes the work of the abductor muscles by retaining corporal heat.
- Recommended for use in rehabilitating processes and in physical exercises.
- Made up of 2 mm. neoprene.

## **Sport Vest**

Code - 872





### **MEASUREMENT PANEL FOR ELASTIC TRUSSES AND OTHER PRODUCTS**

### Code 601\*- 602

Waist Circumference	65-80	80-100	100-125
Corset Nr.	S (1)	M (2)	L (3)

### **Code 701**

Waist Circumference	60-70	70-80	80-90	90-100	100-110	110-125
Suspensory Nr.	5	6	7	8	9	10

### Code 871 - 872

0000011 012						
Waist Circumference	60-70	70-80	80-90	90-100	100-110	110-125
Short Nr.	1	2	3	4	5	6

\* Right - Left



# NE US Unisize / Unisex







Compression Adjustable Products

**NEUS** products are one size. Thanks to its special design and fabric type, it is adjustable on the body.

**NEUS** porducts are made from neoprene, and help to prevent and treat injuries such as pulled muscles and strained tendons.

UNIFORM COMPRESSION: Neoprene provides compression supporting to damaged tissue. Unlike other elastic fabrics it stretches in all directions, but does not overstretch.

**HEAT RETENTION:** Neoprene is a superb insulation which holds radiated body heat onto the affected area stimulating greater blood circulation and aiding a swift recovery.

**COMFORT:** Nexus products are specially designed in order to fit on to the body snugly.




# Unisize / Unisex NE US

# Code - 881 • Helps to relieve pain of injured wrists through pressure.

# **Universal Size**

### **Compression Adjustable Products**



Code - 882

Provides comfortable • compression and supports for weak or stressed wrists.








**Epicondylitis Bandage** 

 Prevention and treatment of epicondylitis while at work or during sport activities.

**Wrist Band** 





Code - 880

Used for holding the toe • static in one position.





## Code - 879

**Wrist Brace Splint** 

- Wrist sprains, tendonitis, tendon synovitis and processes that requires a strong holding of the wrist.
- Traumatic sequelas.
- After orthopedic or surgical treatment.

# NE US Unisize / Unisex





## Code - 884

## **Knee Support**

 Useful in therapy requiring heat and compression.



### **Compression Adjustable Products**

- Spiral stays are put in the lateral grooves intending to provide the best support.
- The front two grooves is used for the thin knees, and the back two grooves is used for the thick knees.





### Knee Support with Flexible Stays

Arthritis, sport injuries, ligament injuries. and for incidents which needs medial-lateral and patellar support.





# Code - 885

# **Open Knee Support**

• This support offers full circumference patellar control for chondromalacia, patellar tracking abnormalities and patellar tendonitis.



### Code - 886

### **Patella Band**

• Especially during the sport activities can be used as a protector for patella.



# **Hinged Stabilizing Knee Brace**

### Code - 894

- Recommended for use in processes of arthrosis, arthritis, bursitis and tendonitis.
- Two aluminium stays, one at right side and one at left side, supports the knee well.



# Unisize / Unisex NE US



**Compression Adjustable Products** 



# Code - 887

# **Ankle Support**

 Provides compression for ankle sprains and soft tissue.



# Ankle Brace With Flexible Stays

### Code - 897



Provides. compression for ankle sprains and soft tissue problems.

Flexible stays at both • sides gives extra support.



# Calf Support

• Used as a support in physiotherapeutic treatment and supports in sport



# Thigh Support Code - 888

Used as a support in • physiotherapeutic treatment and supports in sport activities.



# NE US Unisize / Unisex







# Code - 889

# Waist Support

• Used as a protector for the muscles on the belly and waist.



# **Practical Waist Support**

Code - 890

- It is used as a support to the waist muscles and spine for light waist pains in daily life.
  - Used as a protector fort he muscles on the waist during sport activities.
- Thanks to its special design pelotte, it makes massage effect, and provides comfort feeling.
  - Compression can be adjusted thanks to its extra belt.

# AWARD WINNING product range.





# New Generation Bandages & Support





**Elbow Support** 



1					
_	Elbow Circumference (cm)	21-24	24-27	27-30	30-33
	Elbow Support Nr.	S (1)	M (2)	L (3)	XL (4)







### Advantages of aloe vera microcapsules

• Regeneration, • Neutral pH, • Therapeutic properties.



### Code - 722 Wrist & Thumb Support



^ <del>`</del>			
Wrist Circumference (cm)	13-16	16-19	19-23
Wrist Support Nr.	S (1)	M (2)	L (3)

• Used for support wrist against soft tisue problems.



Code - 723	Calf Supp

Softsport



1	( ) —				
	Calf Circumference (cm)	32-36	36-40	40-44	44-48
	Calf Support Nr.	S (1)	M (2)	L (3)	XL (4)

•Used for muscle injuries, pulled muscles and tendom problems.



Code - 724 Ankle Support



K	- Al				
	Ankle Circumference (cm)	15-19	19-23	23-28	28-33
	Ankle Support Nr.	S (1)	M (2)	L (3)	XL (4)

•Used for ankle sprains and soft tissue problems.







	<b>)</b>				
_	Knee Circumference (cm)	30-34	34-38	38-43	43-48
	Knee Support Nr.	S (1)	M (2)	L (3)	XL (4)

• Used for knee injuries and soft tissue problems.



# V-TEKS Orthopedic Pillows & Cushions VARITEKS ORTOPEDI SANAYI A.Ş.



Code - 651

VISCO Orthopedic Neck Pillow



Available in 3 Sizes:

Small (S) - 47x30x10-8 cm Medium (M) - 51x37x11-8cm (L) - 60x40x11-8cm Large



Code - 652

VISCO Orthopedic Neck Pillow - Classic



Available in 2 Sizes:

Classic - 55x40x10 cm Classic-Kid Size - 42x30x7 cm





Code - 653

VISCO Orthopedic Neck Pillow - Anatomic



53x38x12-9 cm



# VARITEKS Orthopedic Pillows & Cushions V-TEKS

Code - 655

**Code - 657** 

Visco



**Cylinder Pillow** 

**Back Cushion** 

(With Belt)



Code - 675

Polyurethane

Code - 677

Polyurethane







Kod - 678 Polyurethane



Code - 661 Visco

**Anatomic Seating Cushion** 

Code - 681 Polyurethane



# V-TEKS Orthopedic Pillows & Cushions VARITEKS ORTOPEDI SANAYI A.Ş.



Code - 662 Visco

**Wheelchair Seating Cushion** 



Kod - 659

Code - 660

Seating Donut Cushion "D" Form

Code - 679 Polyurethane

Code - 680

Polyurethane



Seating Donut Cushion "U" Form

### NOTES ABOUT PILLOW and CUSHION RANGE:

- The pillows are produced in two groups: VISCO and Polyurethane.
- Orthopedic neck pillows are produced only from VISCO material. Other pillows and cushions are produced from either VISCO and POLYURETHANE material.
- Inner material of the pillows are 100% VISCO or 100% Polyurethane depending on group.
- Inner cover of the pillows and cushion can be washed thanks to the special design.
- Outer cover of the pillows and cushion are made from
- •BAMBOO fabric. The outer part of the cover touching to the skin is 100% BAMBOO.

### **ADVANTAGES of BAMBOO FABRIC:**

• Always keeps the skin dry, for it has 3-4 times more air permeability than cotton.

• Provides feeling of extra softness.

• Naturally antibacterial.











- Strict immobilization after traumatism of cervical rachis.
- Post-surgical immobilization.
- Serious cervical sprains.
- Treatments, which require strict cervical immobilization.
- Thanks to the chin and occipital supports, the flexion and stretch movements of the neck are mechanically restricted and to a lesser degree the lateral and rotational movements It is
- Made up of 12 mm plastazote.







# VARITEKS' kids



- Conservative treatment for elbow dislocations.
- Rehabilitation purpose for elbow contracture.
- Orthotic treatment on the medial or lateral epicondyle fractures.
- Postoperative care after prosthesis implantation.
- Conservative treatment and rehabilitation purpose for elbow fracture.
- Conservative treatment and rehabilitation purpose during postoperative period of elbow ligament injuries.
- Keeps the elbow at an angle or in a certain range of angles.



Code - 859 Elbow ROM Orthosis





Wrist Support

Code - 882



 Can be used for both sides, right and left.



Code - 303 Arm Sling (Sponge)



 Wrist sprains, tendonitis tendon synovitis and processes that requires a strong holding of the wrist.

- Traumatic sequelas.
- After orthopedic or surgical treatment.

Wrist Brace Splint

Code - 825



Code - 894 Hinged Stabilizing Knee Brace

# NE US®

Recommended for use in processes of arthrosis, arthritis, bursitis and tendonitis.



Ankle Support Code - 887

Provides compression • for ankle sprains and soft tissue.



# NE US®

Code - 897

# **Ankle Brace With Flexible Stays**



- Provides compression for ankle sprains and soft tissue problems.
- Flexible stays at both sides gives extra support.

Code - 320 Knee Immobilizer



- Postoperative medial ligament repair, immobilization in the acute phase of corrective osteotomy.
- Assisted walking in operated fractures of the patella.
- Pre-rotulian bursitis.
- In general treatment which requires temporary immobilization of the straightened knee.
- Adjustable in width 20 cm, 25 cm, 30 cm, 35 cm available in length according to the leg.

**Open Knee Support** Code - 885



This support offers full circumference patellar control for chondromalacia, patellar tracking abnormalities and patellar tendonitis.







### **Abduction Corset**

- Used for dislocated hips.
- Made up of 12 mm plastazote.



0000 100					
Waist Circumference (cm)	26-32	32-37	37-42	42-47	
Product Nr.	1	2	3	4	ĺ

## Pavlic Bandage

47-53 5


### Code - 119

- Used for treatments of dislocated hips.
- Provides the baby pull his/her legs towards to body, and separating the knees from each other.
  - Adjustable thanks to the easy buckle system.



25-35 36-46 47-57





# Code - 604 Child Umbilicus

- Used for case of hernia umbilicus at babies.
- It is placed on the belly bottom of baby.
- Soft cushion provides the acceptable pressure on the belly of the baby.



Universal Size


## Elastic Truss For Child Code - 603

- Used for hernia at children.
  - Right and left types are available.









Year 1978 STANDARD KORSTEX

products have accomplished the success of being the first orthopaedic products that were exported from Turkey.







**ISO 13485** certificates, which have been gotten and continuosly being renewed for

VARİTEKS ORTOPEDİ SANAYİ A.Ş

& VARİTEKS MEDİKAL SAN. ve TİC. LTD. ŞTİ.

Our **( 6** declaration for all of our products.

