














	VARICOSE STOCKINGS _____ ●	2-9
	CERVICAL COLLARS _____ ●	10-12
	WAIST CORSETS _____ ●	13-16
	ABDOMINAL & WOMAN CORSETS _____ ●	17-20
	LIPOSUCTION CORSETS _____ ●	21-23
	KNEE BRACES _____ ●	24-27
	ARM SLINGS _____ ●	28-29
	WRIST BRACES & SPLINTS _____ ●	30-32
	SHOULDER & ELBOW SUPPORTS _____ ●	33-34
	ANKLE BRACES & FOOT SUPPORTS _____ ●	35-38
	Silicopad SILICONE PRODUCTS _____ ●	39-43
	ELASTIC TRUSSES AND OTHER PRODUCTS _____ ●	44-45
	NEXUS RANGE (UNISIZE & UNISEX AND ADJUSTABLE PRODUCTS) _____ ●	46-50
	Softsport NEW GENERATION BANDAGES & SUPPORT _____ ●	51
	ORTHOPEDIC PILLOWS & CUSHIONS _____ ●	52-54
	KID SIZE PRODUCTS _____ ●	55-59

* All rights of the photographs reserved. The photographs can not be used without permission.

* It must be taken into consideration that the models, colors and the measurements of the products in the catalog may be changed.




VARICOSE STOCKINGS



All of our varicose stockings are produced according to ISO13485 and CE with the state of the art machines. Best raw materials which are produced for medical stockings in the world are used in all of our stockings. Designed for both men and women.

VARICOSE STOCKINGS

Code - 902 Knee High (Open Toe)



High Compression

**30-40 mm.Hg
Class III**

- Constitution arising from constitutional or post-thrombotic venous insufficiency.
- Strong tendency to edema.
- Secondary varicosis.
- White atrophy, dermatosclerosis.
- Following treatment for severe ulcers, if they have caused a relapse.


Code - 901 Thigh High (Open Toe)



Silicone Dots



Code - 900 Metal Stocking Donning Kit



Code - 963 Pantyhose (Open Toe)

NEW



- **Wearing Satin** is put on the foot for open toe before wearing the stocking.
- Stocking can be worn easily.
- It can be removed holding from toe part.

VARICOSE STOCKINGS

Code - 912 Knee High (Open Toe)



Advantage of aloe vera microcapsules

- Regeneration, • Neutral pH, • Therapeutic properties.

Medium Compression

**20-30 mm.Hg
Class II**

- Varicosis marked by a tendency to edema.
- Post-traumatic swelling.
- Following treatment for light ulcers and superficial thrombophebitis.
- Following sclerisin and surgical treatment of varicose veins.
- Severe pre-natal varicosis.

Advantages of Microfiber:

- Microfiber is super soft to wear.
- Gives full breathability.
- Removes moisture as it is produced.



Code - 911 Thigh High (Open Toe)

Code - 913 Pantyhose (Open Toe)



Code - 917 Pregnancy Pantyhose (Open Toe)

NEW



- **Wearing Satin** is put on the foot for open toe before wearing the stocking.
- Stocking can be worn easily.
- It can be removed holding from toe part.



Available in Beige and Black colors.



VARICOSE STOCKINGS



Medium Compression

**20-30 mm.Hg
Class II**

Knee High (Closed Toe)

Code - 905



Code - 904

Thigh High (Closed Toe)



Code - 903

Pantyhose (Closed Toe)



Code - 907

**Pregnancy Pantyhose
(Closed Toe)**



Code - 920

Varicose Socks For Men*



* Available in Black and Dark Grey colors



Available in Beige and Black colors.

VARICOSE STOCKINGS



Light Compression

**15-20 mm. Hg
Class I**

- Symptoms of spider veins. Minor swelling, achy legs.
- Small varicose veins and post operation.

Code- 931 Thigh High (Closed Toe)

Code - 932 Knee High (Closed Toe)

NEW



Code - 933 Pantyhose (Closed Toe)

NEW



Super Soft
Multifilament

Advantages of Multifilament:

- Supersoft to wear.
- Fits well thanks to the elasticity.



Available in Beige and Black colors.

VARICOSE STOCKINGS

Code - 910 Knee High



**NEW
SYSTEM**

Code - 909 Thigh High



Anti-Embolism

15-20 mm.Hg

- Clinically designed to promote venous return for prevention of thrombosis before, during and after surgery and for medical conditions requiring prolonged bed rest and before, during and after child birth.

Code - 908 Thigh High With Belt



Super Soft
Multiflament

Advantages of Multiflament:

- Supersoft to wear.
- Fits well thanks to the elasticity.

VARICOSE STOCKINGS

Code - 803 Elastic Knee Brace



- Used for all knee pains, injured knees and pulled muscles of knee.

Elastic Tendon Brace Code - 802

- Used for muscle injuries on the lower leg, for pulled muscles and circulation difficulties.
- Also for those who can not wear the ankle part of varicose stocking it serves as such.



Code - 801 Elastic Ankle Brace

- Used for all pains around the ankles, looseness around the ankles.



NEW



Medium Compression

**20-30 mm.Hg
Class II**

- Postoperative and posttraumatic lymphedema treatment.
- Especially for treatment of edemas after mastectomy operations.

NEW




Code - 951 Lymphedema Arm Sleeve (With Belt)

Lymphedema Arm Sleeve Code - 952


VARICOSE STOCKINGS

MEASUREMENT PANEL FOR VARICOSE STOCKINGS


Code 901 - 902 - 963 - 903* - 904* - 905* - 907* - 911* - 912* - 913* - 917* - 802

	Calf Circumference (cm) cC	30-32	32-34	34-36	36-38	38-40	40-42	42-45
	Stocking Nr.	1	2	3	4	5	6	7


Code 931* - 932* - 933* - 920**


	Calf Circumference (cm) cC	29-32	32-35	35-38	38-41	41-44
	Stocking Nr.	S (1)	M (2)	L (3)	XL (4)	XXL (5)

Code 801

	Ankle Circumference (cm) cB	15-17	17-19	19-21	21-23	23-25	25-27	27-30
	Ankle Brace Nr.	1	2	3	4	5	6	7

Code 803

	Knee Circumference (cm) cE	28-31	31-34	34-37	37-40	40-43	43-46	46-50
	Knee Brace Nr.	1	2	3	4	5	6	7



* Available Beige and Black colors.

** Available Black and Dark Grey colors.

Anti Embolism Measurement Panel


Product Nr.	S (1)	M (2)	L (3)	XL (4)	XXL (5)
Calf (cC) Circumference (cm) 908-909-910	28-32	32-36	36-40	40-44	44-49
Leg (IG) Length (cm) 908-909	Regular	up to 80	up to 80	up to 80	up to 80
	Long	over 80	over 80	over 80	over 80
Knee (IE) Length (cm) 910	Regular	up to 40	up to 42	up to 45	up to 48
	Long	over 40	over 42	over 45	over 48

Note: Code 908 is available in right and left sides.

On the stockings;

Colors on the toe show the sizes: **Green S(1)**, **Yellow M(2)**, **Red L(3)**, **Blue XL(4)**, **Grey XXL (5)**


Colors on the upper side show the length. **Yellow Regular**, **Blue Long**



Lymphedema Arm Sleeve Measurement Panel

Code 951 - 952

	1	2	3	4	5	6
Circumference Ref. cC (cm)	14,5-16	16-17,5	17,5-19	19-20,5	20,5-22,5	22,5-24,5
cC1 (cm)	17-18,5	18,5-20	20-22	21,5-24	23,5-26,5	26-29
cE (cm)	24,5-28	26,5-29	28-31	29,5-32,5	31-34	33-37
cG (cm)	26-31	29-34	31-36	33-38	35-41	37-43



CERVICAL COLLARS



- Used for broken necks, burns and after neck operations.
- Long velcro tab for universal application and adjustment.

Code - 201 Support Collar



- Used for broken necks, burns and after neck operations
- Long velcro tab for universal application and adjustment

Chin Support Collar Code - 202




Code - 203 Soft Collar

- Comfortable foam (45 density) covered with stockinet, which is made from **BAMBU** yarn absorbing bothersome perspiration. Used for slight cervical problems.



Nelson Collar Code - 204

- Used for soft tissue damages and osteoarthritis
- Practical useable
- It is made up of 12 mm. plastazote



NEW



- Used for soft tissue damages and slight cervical problems.
- It is made from comfortable foam fabric which is in optimum strength.
- Height adjustment can be done thanks to the velcro which comes under chin.

Universal Boy

Nr. cm:	(S) 25-30	(M) 30-35	(L) 35-40	(XL) 40-46
---------	-----------	-----------	-----------	------------

UNI-COLLAR Soft

Code - 208

CERVICAL COLLARS



Code - 206 Nelson Collar With Chin Support

- Used for broken necks, burns and after neck operations.
- Keeps the neck in extension position.
- It is made up of 12 mm. plastazote



Soft Nelson Collar Code - 207

- Used for soft tissue damages and osteoarthritis
- Practical useable
- It is made from 12 mm. plastazote (inside)
- Covered by soft sponge fabric.



Code - 205 Philadelphia Collar

- Strict immobilization after traumatism of cervical rachis.
- Post-surgical immobilization.
- Serious cervical sprains.
- Treatments, which require strict cervical immobilization.
- Thanks to the chin and occipital supports, the flexion and stretch movements of the neck are mechanically restricted and to a lesser degree the lateral and rotational movements.
- It is made up of 12 mm. plastazote.



Philadelphia Tracheostomy Collar Code - 210

- Used for all the events for which philadelphia collar is used, and for tracheostomy events.

CERVICAL COLLARS



Code - 260 **Nose Guard Mask**

- It is used in order to protect the nose after rhinoplasty operations.
- It is a protector during some sports where the helmet is not used.



Minerva Collar **Code - 250**

- Designed to keep patient's head in an immobilized position which requires fast and secure immobilization
- Dorsal section clearly allows the patient to lie flat without discomfort while the head remain in prescribed position.

- If chin to umbilicus under 40 cm - SMALL
- If chin to umbilicus over 40 cm - BIG



MEASUREMENT PANEL FOR CERVICAL COLLARS

Code 201 - 202 - 203

Neck Circumference (cm)	22 - 27	27 - 32	32 - 38	38 - 44	44 - 50
Collar Nr.	S (1)	M (2)	L (3)	XL (4)	XXL (5)

Code 204 - 205 - 206 - 207 - 210

Neck Circumference (cm)	25 - 30	31 - 38	39 - 46
Collar Nr.	S (1)	M (2)	L (3)

Note: 203 - 204 - 205- 206 are available in kid sizes.



WAIST CORSETS



Code - 140

Elastic Dorsolomber Corset
(Height Adjustable)

- For disk hernia, displaced disks, misshaped disks, scoliosis, kyphosis broken spine and during recovery from these ailments.
- Height adjustable thanks to the two pieced system.

- For disk hernia, displaced disks, misshaped disks, scoliosis, kyphosis, broken spine and during recovery from these ailments
- Metal parts of the corset are made up of hard stainless steel, which prevents deformation



Dorsolomber Corset

Code - 110



Code - 130

Elastic Dorsolomber Corset

- For disk hernia, displaced disks, misshaped disks, scoliosis, kyphosis broken spine and during recovery from these ailments.



TX Dorsolomber Corset

Code - 162

- For disk hernia, displaced disks, misshaped disks, scoliosis, kyphosis broken spine and during recovery from these ailments.
- Fits to the body well and full support thanks to the elastic laces and belt.
- Gives breathability thanks to the special design fabric.

WAIST CORSETS



- For disk hernia, lumbago and sciatica pains.
- Supports the back, and prevents body weight to fall on the small of the back.
- Metal parts of the corset are made up of hard stainless steel, which prevents deformation.



- For disk hernia, lumbago and sciatica pains.
- Supporting the back, prevents body weight to fall on the small of the back.
- 31 cm in height.
- Fits to the body well, and full support thanks to the elastic laces and belt.
- Gives breathability thanks to the special design fabric.

Code - 111 Lumbosacral Corset

TX Lumbosacral Corset Code - 161



- For disk hernia, lumbago and sciatica pains.
- Supporting the back, prevents body weight to fall on the small of the back.



Code - 135 Elastic Lumbostad Corset



- 28 cm in height.



- For lumbago, sciatica, disk hernia and displaced disk, which are not severe.

- 20 cm in height.

Elastic Light Lumbostad Corset Code - 145

WAIST CORSETS



- For lumbago, sciatica, disk hernia and displaced disk which are not severe.
 - Made up of elastic cloth.
 - Easy to use thanks to the elastic fabric and hook and loops.
- Another model with rigid steels is available. (KOD115-C)

Code - 115 Lumbostad Corset


- Suitable for lumbar and rheumatic pains and for prevention when the lumbar area is subject to continued stress.
- Made up of 3 mm neoprene.



Neoprene Lumbostad Corset Code - 821



NEW



S	M	L	XL
65-75 (cm)	75-90 (cm)	90-105 (cm)	105-125 (cm)

- Size adjustment is made according to the above table.
- Size adjustment is performed by moving the right and left sides of the corset among the coloured area.
- Right and left sides both must be stabilized on the same coloured part.
- Right and left sides can be replaced each other depending on patients usage direction.



Universal Size

Code - 181 Sacro-Plus Back Brace PRO



NEW



Sacro-Plus Back Brace

Code - 185

- Helps treatment by supporting lumbar and sacral area at the events of disk hernia, lumbago and sciatica pains.
- The back support mechanism, provides optimum stabilization to the disks, and also helps the corset to fit well on the waist according to vertebra anatomy.
- By supporting the waist, it reduces the effect of the body weight falling on the small of the back.

WAIST CORSETS



MEASUREMENT PANEL FOR WAIST CORSETS

Code 110 - 111

Waist Circumference(cm)	60-70	70-80	80-90	90-100	100-110	110-120	120-130
Corset Nr.	1	2	3	4	5	6	7

Code 130 - 131 - 135 - 138 - 140 - 145 - 821

Waist Circumference(cm)	60-70	70-85	85-100	100-115	115-130
Corset Nr.	S (1)	M (2)	L (3)	XL (4)	XXL (5)

Code 118

Waist Circumference(cm)	65-75	75-85	85-95	95-105	105-120
Corset Nr.	S (1)	M (2)	L (3)	XL (4)	XXL (5)

Code 115 - 115C

Waist Circumference(cm)	70-75	75-80	80-85	85-90	90-95	95-100	100-105	105-110	110-120	120-130
Corset Nr.	65	70	75	80	85	90	95	100	105	110

Code 161 - 162

Waist Circumference(cm)	60-75	75-90	90-110	110-130
Corset Nr.	S (1)	M (2)	L (3)	XL (4)

Note: Code 118, 130 and 138 are available in kid size.



Code - 118 Clavicle Bandage

- Clavicle band gives accurate immobilization and prevents dislocation.
- It provides a warning for incorrect posture.



Posturex (Clavicle) Corset Code - 138

- Prevention of juvenile khyphosis.
- Rehabilitating treatment of dorsum curve or weighted back.
- Ideal for keeping the back straight, avoiding incorrect posture.



ABDOMINAL & WOMAN CORSETS



• 28 cm in height.



Code - 132

Postoperative Abdominal Corset (28 cm.)

- Thoracic and abdominal injuries.
- Post-operative recovery periods.
- Post-natal periods, abdominal hernias.
- Abdominal muscular traumatism.



20 cm in height. •



Postoperative Abdominal Corset (20 cm.)

Code - 133



NEW

Code - 136 Colostomy Corset

- It is used in order to reduce the risk of parastomal hernia after colostomy operations.
- Specially designed to support the colostomy bag and abdominal area.



NEW

Code - 142 Abdolumbar Korset

- Thoracic and abdominal injuries.
- Post-operative recovery periods.
- Post-natal periods, abdominal hernias.
- Abdominal muscular traumatism.
- It has four pieces of plastic stays in order to support the waist and to prevent the corset to curl back.
- 28 cm in height.



Elastic Abdominal Hernia Corset Code - 134

- Used for abdominal hernia in order to prevent it to grow.
- 20cm in height.

ABDOMINAL & WOMAN CORSETS



JUDO[®]

- Used in order to keep the waist, abdomen and kidneys warm.
- For spasms of muscles.
- For pains due to rheumatism and sciatica.
- Blend of angora and wool.

Code - 114 Angora Waist Corset



JUDO[®]

- To keep the waist, abdomen and kidneys warm
- For spasms of muscles.
- For pains due to rheumatism and sciatica.

Waist Corset Code - 112



Code - 101 Hooked Waist Corset

- For back pains and to shape the body.
- Abdominal muscular traumatism.



Panty Corset Code - 113

- For rheumatism, sciatica, lumbago, kidney infection.
- For infection of intestines and gynecological diseases.

JUDO[®]

ABDOMINAL & WOMAN CORSETS



Code - 137 Rib Belt Corset

- Used after all chest operations or any complications involving the chest.

- Post pregnancy deformations of the abdomen, deformations, after weight loss, abdominal hernias.
- Elasticity and the special design provides comfort.
- Made from elasthan.



Corset Type Skandal Code - 105



- Post pregnancy deformations of the abdomen, weight loss, uterus misplacements and hernias.
- Made from elasthan.

Code - 106 Corset Type Panty

- Helps to lose weight, and helps dissolve fat tissue around hips and abdomen.
- It is also used after gynecological operations and liposuction operations.
- Made from elasthan.



Sanitary Panty Corset Code - 107

ABDOMINAL & WOMAN CORSETS



- This corset has two functions:
 - 1 Reduces probability of miscarriage before the child birth.
 - 2 After the child birth, it is reversed, and used, like in the small photograph, in order to prevent sagging and deformation.
- Removes waist pains. It helps the fetus to rest comfortably in the uterus.

Kod - 139 Elastic Pregnancy Corset



- Used before the child birth.
- Prevents sagging and cracking events.
- It helps the fetus to rest comfortably in the uterus.
- Reduces the risk of miscarriage.
- Adjustable with the laces during the period of pregnancy.



Pregnancy Corset Code - 109

MEASUREMENT PANEL FOR ABDOMINAL & WOMAN CORSETS

Code 101

Waist Circumference(cm)	65-70	70-75	75-80	80-85	85-90	90-95	95-100	100-105	105-110	110-120
Corset Nr.	60	65	70	75	80	85	90	95	100	105

Code 105 - 106 - 107

Waist Circumference(cm)	60-65	65-75	75-85	85-95	95-105	105-120
Corset Nr.	S (1)	M (2)	L (3)	XL (4)	XXL (5)	XXXL (6)

Code 109 - 139

Body Size Nr.	36-38	40-42	44-46	48-50	52-56
Corset Nr.	S (1)	M (2)	L (3)	XL (4)	XXL (5)

Code 112 - 113 - 114 - 132 - 133 - 134 - 136 - 142

Waist Circumference(cm)	60-70	70-85	85-100	100-115	115-135
Corset Nr.	S (1)	M (2)	L (3)	XL (4)	XXL (5)

Code 137

Chest Circumference(cm)	65-75	75-85	85-95	95-105	105-120
Corset Nr.	S (1)	M (2)	L (3)	XL (4)	XXL (5)



LIPOSUCTION CORSETS



NEW

Code - 235 Liposuction Chin - Neck Bandage

NEW MATERIAL



NEW

Liposuction Corset
(Chest - Arm - Back)

Code - 230

Liposuction Corsets

- Liposuction corset is a modern compressive outfit, ideal wear for the post-operative care following plastic surgery, such as liposuction around the waist, hip, abdomen, calf, arm and neck.
- Provides the skin to breath thanks to its special fabric.
- Provides the body micro-massage effect thanks to the four-way stretch fabric.



Code - 129 Liposuction Arm Corset



Liposuction Panty Corset Code - 124



Available in Beige and Black colors.

LIPOSUCTION CORSETS



Code - 121

**Liposuction Corset
(Under Breast - Below Knee)**




Code - 125

**Liposuction Corset
(Under Breast - Ankle)**




**Liposuction Corset
(Under Breast - Above Knee)**

Code - 122



Code - 123

**Liposuction Corset
(Waist - Below Knee)**



**Liposuction Corset
(Waist - Ankle)**

Code - 126

Available in Beige and Black colors.

LIPOSUCTION CORSETS

NEW



- Used by men after gynecomastia operations.
- Covers the operated area well.
- It is easy to apply the corset thanks to the hook and zipper system.



Easy buckle system provides comfort.



Available in Black and Beige colors.




Code - 240

Liposuction Men's Gynecomastia Vest

MEASUREMENT PANEL FOR LIPOSUCTION CORSETS

Code 121 - 122 - 123 - 124 - 125 - 126

 Waist Circumference(cm)	60-65	65-75	75-85	85-95	95-105	105-120
Corset Nr.	S (1)	M (2)	L (3)	XL (4)	XXL (5)	XXXL (6)


Code 129

 Upper Arm Circumference (cm)	23-26	27-29	29-33	33-37	37-41
Corset Nr.	S (1)	M (2)	L (3)	XL (4)	XXL (5)

Code 230

 Under Breast Circumference (cm)	67-72	73-79	80-88	89-97	98-106	107-115
Corset Nr.	S (1)	M (2)	L (3)	XL (4)	XXL (5)	XXXL (6)

Code 240

 Chest Circumference (cm)	92-96	96-102	102-110	110-118	118-126	126-135
Corset Nr.	S (1)	M (2)	L (3)	XL (4)	XXL (5)	XXXL (6)

Code 235

 Neck Circumference (cm)	24-35	34-46
Corset Nr.	S	L